

## Continuous Processing

From process development to customized production, continuous separation and work-up

### Evonik has more than 50 years of experience in developing and applying continuous processes:

- Established processes with micro reactors, tube reactors, loop reactors, trickle-bed reactors, and many more.
- In-house reactor optimization / development.
- Continuous separation and work-up technologies: Distillation, extraction, crystallization, membranes.
- Simulation and design capabilities at our Hanau, Germany site.

### Pilot plant

A flexible modular continuous processing pilot plant at our site in Hanau, Germany.

- Main reaction types: Hydrogenation, high and low T, hazardous gases
- Upgrade to cGMP possible
- OEB 3


### Large-scale production

- More than 60 % of Evonik's products are produced using continuous processes.
- In-house process technology and engineering unit (> 500 engineers and chemists) for complex projects.
- Track record of complex and large investment projects (also in exclusive customer projects).

### Big engineering network

- Experienced engineering capabilities for special topics
- Cooperation with experts:
  - Process Analytic Technology (PAT) group
  - Simulation and additive manufacturing (3D-printed reactors)

### Example 1: Continuous hydroxyalkylation in a tube reactor



### Significant improvements by switching from batch reaction to tube reactor:

- Yield of batch process was improved by 50 %
- Consistently higher quality of final product, fewer impurities
- More than 10 campaigns have been run so far at Evonik

## Process development capabilities

### Modular Small-Scale Flow System:

- Different mixers/reactors (tube, micro mixer, static mixer, T-mixer, dynamic magnetic driven mixer)
- -70 to +200 °C
- Residence times from 3 seconds to 30 minutes
- Flow rate: 6 to 1200 ml/h
- Two consecutive reaction units with independent temperature control

### Loop/Tube Reactor:


- -40 to +170 °C
- Residence times from 3 to 30 minutes
- Flow rate: 50 to 500 ml/h
- Two consecutive reaction units with independent temperature control

### Continuously Stirred Reactor Cascades:


- Automated parallel equipment: 30 – 100 ml
- 1 – 2 l double jacketed vessels
- Cascades of up to 8 reactors
- -10 to +140 °C
- Residence time: 10 – 120 min/stage

### Process analytical technologies:

- Collaboration with Process Analytical Technology (PAT) group
- PAT technology
  - FTIR
  - RAMAN
  - Automated sampling



## Example 2: Continuously stirred tank reaction to produce an API intermediate (600 kg scale)



- Grignard run as a continuously stirred tank reaction in 100 l reactor
- Manufactured several tons of Grignard solution

This information and all further technical advice are based on our present knowledge and experience. However, it implies no liability or other legal responsibility on our part, including with regard to existing third-party intellectual property rights, especially patent rights. In particular, no warranty, whether expressed or implied, or guarantee of product properties in the legal sense is intended or implied. We reserve the right to make any changes according to technological progress or further developments. The customer is not released from the obligation to conduct careful inspection and testing of incoming goods. Performance of the product described herein should be verified by testing, which should be carried out only by qualified experts in the sole responsibility of a customer. Reference to trade names used by other companies is neither a recommendation, nor does it imply that similar products could not be used.

**Evonik Operations GmbH**  
Health Care Business Line

healthcare@evonik.com  
www.evonik.com/healthcare