Technological Capabilities

Organometallic Chemistry and Cryogenic Capacities

Decades of experience in boronic acids, highly selective reductions, and Li/Mg/B chemistry

Setup

- The largest cGMP capacity for cryogenic reactions:
- > 200 m³ cryogenic reactor volume
- Reactions down to -80 °C
- 500 l to 8 m³ reactors
- Glass lined and Hastelloy
- Up to 10 bar
- OEL down to 0.1 μg/m³
- · Batch mode and continuous processing
- Highly flexible setup with more than 50 reactors with cryogenic capacities
- · Organometallic reactions at three sites



Reaction types

- Decades of experience in organometallic and cryogenic chemistry including
- · Grignard reactions
- · Lithiations: LDA, BuLi, HexLi, etc.
- Reductions: LAH, DIBAH, BH₃, 9-BBN, etc.
- Aromatic/heterocyclic boronic acids
- Sensitive reagents (ZnEt2 etc.)
- Transmetallations
- · Highly exothermic reactions
- · Chemistry in Water

Special technologies

- Grignard in continuously stirred tank reactor cascades
- Loop reactors
- Large cryogenic high potent capacities
- Precious Metal Recovery at Tippecanoe
- Sustainability Aspects

Example 1: Boronic acid formation in a loop reactor

$$\begin{array}{c} X \\ R \end{array} \qquad \begin{array}{c} LDA \\ (^{i}PrO)_{3}B \\ \hline THF -10 \text{ to } 0^{\circ}C \end{array} \qquad \begin{array}{c} O^{i}Pr \\ I \\ B \\ \end{array} \qquad X \\ R \end{array} \qquad \begin{array}{c} OH \\ HO^{-B} \\ X \\ R \end{array}$$

- Cost Improvements via Precious Metal Recovery at Tippecanoe for large scale productions
- Membrane Filtration Technology enables Recovery of precious Metals from organic waste streams
- Sustainable Technology with Significant Cost Improvement



Example 2: Rearrangement followed by Negishi coupling

- Handling of Grignard and organic Zinc compounds (implemented at ton scale)
- · Highest purity requirements
- Involves continuous step
- More than 10 years experience in commercial production

Example 3: Continuously stirred tank reaction to produce an API intermediate (600 kg scale)

- Grignard run as continuously stirred tank reaction in 1001 reactor
- Manufactured several t of Grignard solution

This information and all further technical advice are based on our present knowledge and experience. However, it implies no liability or other legal responsibility on our part, including with regard to existing third-party intellectual property rights, especially patent rights. In particular, no warranty, whether expressed or implied, or guarantee of product properties in the legal sense is intended or implied. We reserve the right to make any changes according to technological progress or further developments. The customer is not released from the obligation to conduct careful inspection and testing of incoming goods. Performance of the product described herein should be verified by testing, which should be carried out only by qualified experts in the sole responsibility of a customer. Reference to trade names used by other companies is neither a recommendation, nor does it imply that similar products could not be used.

Evonik Operations GmbH Health Care Business Line

healthcare@evonik.com www.evonik.com/healthcare

