

## Transition-Metal Catalysis

Your experienced partner with a full service offering

### > 30 years of experience in developing, tailoring & commercializing catalysts for organic reactions

Our in-house library of homogenous catalysts comprises:

- Approx. 1200 ligands
- Approx. 250 metal precursors
- All important ligand classes.

Catalyst/product separation by nanofiltration for catalyst recycling (Evonik membrane technology).

### > 20 years of experience in catalytic process development

#### A complete portfolio of catalyst technology

Evonik is a leading supplier of heterogeneous catalysts

- Precious metal powder catalysts (> 100 commercial products)
- Custom-made precious metal powder catalysts
- Raney-type catalysts (> 10 commercial products)
- Continuous process catalysts

### Inhouse catalyst production


Evonik offers the advantages of a combination of in-house catalyst development, production, as well as the development and scale up of the final organic process.

#### Your advantage:


- Superior technical solutions
- Quicker timelines
- Secure supply chain
- Seamless catalyst recovery

### Broad range of covered reaction types

- (Enantioselective) hydrogenation
- (Enantioselective) reductive amination
- Hydroformylation
- C-X coupling
- C-C coupling
- Metathesis
- Oxidation
- Dehydrogenation
- and many more



**Route 1:** Hydrogenation of benzaldehyde at < 6 bar. Addition of metal salts allowed for quantitative conversion at 6 bar instead of 20 bar.



**Route 2:** Selective hydrogenation of phenols. Doping of Pd/C catalysts allowed a high conversion of phenol to cyclohexanon with minimal formation of cyclohexanol.

### From first lead detection to a scaleable process – our capacities around the world:

- Globally cGMP production of APIs in five different sites
- Up to 25 bar
- 80 °C – 200 °C
- 100 – 16.000 l
- Glass lined or Hastelloy
- Highly potent compounds



### From first lead detection to a scaleable process

- First lead detection (HTS): 96 well plates; up to 100 °C, 100 bar
- Catalyst optimization (HTS): Four reactor batteries (4 x 8 reactors/2 ml); up to 140 °C, 150 bar
- Process development: Various reactors (30 ml – 300 ml): Up to 140 °C, 150 bar
- Scale up: Various autoclaves: Monel-, stainless steel- or glass-made (1000–2500 ml): Up to 60 bar, 150 °C



This information and all further technical advice are based on our present knowledge and experience. However, it implies no liability or other legal responsibility on our part, including with regard to existing third-party intellectual property rights, especially patent rights. In particular, no warranty, whether expressed or implied, or guarantee of product properties in the legal sense is intended or implied. We reserve the right to make any changes according to technological progress or further developments. The customer is not released from the obligation to conduct careful inspection and testing of incoming goods. Performance of the product described herein should be verified by testing, which should be carried out only by qualified experts in the sole responsibility of a customer. Reference to trade names used by other companies is neither a recommendation, nor does it imply that similar products could not be used.

**Evonik Operations GmbH**  
Health Care Business Line

healthcare@evonik.com  
www.evonik.com/healthcare