

**Outstanding strength. Up
to 6x faster degradation.**


For bioresorbable implant applications

Resomer®

PLA-PEG copolymers



Bioresorbable, biocompatible PLA-PEG triblock copolymers




RESOMER® PEG is a commercial platform of triblock copolymers comprised of hydrophobic polylactide (PLA) blocks and hydrophilic polyethylene glycol (PEG) blocks (PLA-PEG-PLA).

The mechanical properties you need. The rapid degradation you want.

Each PLA-PEG copolymer in the standard range replicates the mechanical strength of its equivalent standard polymer but provides degradation rates up to six times faster. With multiple customization options available upon request, the

platform delivers the versatility to enhance performance across a range of applications including wound healing and pediatric implants. Each PLA-PEG copolymer also leverages more than 30 years of safety data for RESOMER®.



Easy to process, with consistent, high-quality performance

All RESOMER® PEG copolymers are purified for low residual monomer content and are compatible with all

relevant technologies including extrusion, compression and injection molding.

NAME	CONTENT		MECHANICAL PROPERTIES		PHYSICAL PROPERTIES	
	MONOMER	PEG	TENSILE STRENGTH (MPa)	YOUNG'S MODULUS (MPa)	Tm °C	Tg °C
RESOMER® LP t 16	100 % L-lactide	1 %	70–80	3500–4500	177	59
RESOMER® LP t 46	100 % L-lactide	4 %	65–75	2500–3500	172	53
RESOMER® LRP t 7016	70 % L-lactide / 30 % D, L-lactide	1 %	60–70	3500–4500	Amorphous	53
RESOMER® LRP t 7046	70 % L-lactide / 30 % D, L-lactide	4 %	60–70	3000–4000	Not applicable	50

This information and all further technical advice are based on our present knowledge and experience. However, it implies no liability or other legal responsibility on our part, including with regard to existing third party intellectual property rights, especially patent rights. In particular, no warranty, whether express or implied, or guarantee of product properties in the legal sense is intended or implied. We reserve the right to make any changes according to technological progress or further developments. The customer is not released from the obligation to conduct careful inspection and testing of incoming goods. Performance of the product described herein should be verified by testing, which should be carried out only by qualified experts in the sole responsibility of a customer. Reference to trade names used by other companies is neither a recommendation, nor does it imply that similar products could not be used.

Evonik Nutrition & Care GmbH
Health Care Business Line
Pharma Polymers & Services

resomer@evonik.com
www.evonik.com/resomer